QIAcube

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4th Edition

Sample & Assay Technologies

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
3. Enzyme tube (1.5ml micro tube/2.0ml Sample tube CB)

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
General Information (May 2008)

Application	DNA
Kit	QIAamp® DNA Mini Kit (50) plus QIAamp DNA Accessory Set A, cat. no. 1048145 or QIAamp DNA Mini Kit (250) plus QIAamp DNA Accessory Set B, cat. no. 1048146
Sample material	Bacteria (Gram+) or yeast
Short protocol name	Enzymatic lysis
Version	1
Full protocol name	Purification of bacterial or yeast DNA with enzymatic lysis V1
Editable parameters	Incubation time and temperature of enzymatic lysis (default 30 min at 37°C) and elution step 2: 50–250 μ l in increments of 25 μ l; default 150 μ l (see “Comments” on the next page).
Required QIAcube® software versions	Firmware version FIW-50-001-J_FW_MB.hex and PLC program version FIW-50-002-G_PLG_MB.prs or higher; available from the QIAcube Web Portal

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
* Sarstedt, Micro tube 1.5 ml, Safety Cap (see www.sarstedt.com).

Volume of reagent required for the indicated number of samples (µl)			
Number of samples	A	B	C
2	63		
3	85		
4	106		
5	128		
6	150		
7	171		
8	192		
9	214		
10	236		
12	279		

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Shaker

Material	Bacterial or yeast cells resuspended in 180 µl enzyme solution (see “Comments” on the next page)
Vessel	2 ml safe-lock microcentrifuge tube*
Adapter	Shaker adapter for 2 ml microcentrifuge tubes (marked with “2”)

Comments



Things to do before starting

Use a pipet to transfer the required volume of Buffer AL to the reagent bottle. Avoid generating large air bubbles or foam as this may lead to a liquid-detection error during the load check.

Note: Lysis of Gram-positive bacteria or yeast is carried out enzymatically. Before starting the protocol, pellet the bacterial cells (maximum 2×10^9 cells) or yeast cells (maximum 5×10^7 cells), resuspend in 180 µl enzyme solution (not provided in the QIAamp DNA Mini Kit), and start the protocol. The temperature and length of the

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QIAcube
Protocol Sheet

General Information (October 2007)

Application	① RNA
Kit	② RNeasy® Mini
Sample material	③ Animal cells
Short protocol name	④ QIAshredder
Version	3
Full protocol name	Purification of total RNA from $5 \times 10^6 - 1 \times 10^7$ animal cells (with QIAshredder homogenization)
Editable parameters	Elution volume: 30–100 µl in increments of 10 µl; default 50 µl

Shaker

Material	$5 \times 10^6 - 1 \times 10^7$ cells disrupted in 600 µl Buffer RLT (thaw cells before starting)
Vessel	2 ml safe-lock microcentrifuge tube*
Adapter	Shaker adapter for 2 ml microcentrifuge tubes (marked with "2")

* Eppendorf®, Safe-Lock Tube 2.0 ml (see www.eppendorf.com).

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

LoadingChart

Appendix

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Products	Contents	Cat. no.
Accessories		
Starter Pack, QIAcube	Pack includes: reagent bottle racks (3); rack labeling strips (8); 200 ul filter-tips (1024); 1000 ul filter-tips (1024); 1000 ul filter-tips, wide-bore (1024); 30 ml reagent bottles (18); rotor adapters (240); rotor adapter holder	990395
Filter-Tips, 1000 ul (1024)	Sterile, Disposable Filter-Tips, racked; (8 x 128)	990352
Filter-Tips, 1000 ul, wide-bore (1024)	Sterile, Disposable Filter-Tips, wide-bore, racked; (8 x 128)	990452
Filter-Tips, 200 ul (1024)	Sterile, Disposable Filter-Tips, racked; (8 x 128)	990332
Rotor Adapters (10 x 24)	For 240 preps: 240 Disposable Rotor Adapters; for use with the QIAcube	990394
Sample tubes RB	QIAcube 2.0ml tube for sample; (1000/pack)	990381
Sample tubes CB	QIAcube 2.0ml tube for enzyme; (1000/pack)	990382
Reagent bottle	QIAcube reagent bottles; 6 x 30ml	990393