
Product Catalog
Instruments

Affinity. Stability.
BIOPHYSICS in a new dimension

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Enjoy the benefits of MST :

Optimize assays quickly :

- ▶ Judge and improve sample quality immediately

Measure previously unmeasurable targets :

- ▶ Work with very small amounts and sensitive samples

Benefit from close-to-native conditions :

- ▶ Analyze in all buffers and biofluids (cell lysate, serum)
– immobilization-free

Do your research efficiently :

- ▶ Enjoy perfect ease-of-use, purification-free measurements and get rid of maintenance downtime

Work flexibly :


- ▶ K_d s for all molecular weights from ions to ribosomes and for pM to mM binding affinities

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
Monolith[®] NT.115 Series

Product Information

Monolith[®] Instruments
for MicroScale Thermophoresis

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- ▶ Immobilization-free affinity determination from 1 pM to mM
- ▶ Broad application range
- ▶ Buffer independency : including serum or cell lysate
- ▶ Purification free : fluorescent fusion proteins

Cat #	Instrument	Channel 1	Channel 2	Affinity Range (K_d)	Sample Consumption (per K_d)
G006	NT.115 ^{Pico}	Pico – RED	–	1 pM to mM	120 pg ¹
G007	NT.115 ^{Pico} Blue/Green	Nano – BLUE	Nano – GREEN	1 nM to mM	120 ng ²
G008	NT.115 Blue/Red	Nano – BLUE	Nano – RED	1 nM to mM	120 ng ²
G009	NT.115 Green/Red	Nano – GREEN	Nano – RED	1 nM to mM	120 ng ²

¹ calculated for a standard protein of 50 kDa, 12 data points per K_d and 10 pM fluorescently labeled protein.

² calculated for a standard protein of 50 kDa, 12 data points per K_d and 10 nM fluorescently labeled protein.


Technical Details

Monolith Instruments NT,115Series	NT,115	NT,115 ^{Pico}
Samples per run	16 samples	16 samples
Fluorescence channels per instrument	2 (BLUE, RED or GREEN)	1 (RED)
Affinity range	1 nM to mM	1 pM to mM
Labeling required	Yes	Yes
Concentration of fluorescent molecule	$10^{-9} - 10^{-3}$ M	$10^{-11} - 10^{-3}$ M
Range of accessible interactions	■ ■ ■ ■ ■	■ ■ ■ ■ ■
Biophysical parameters	Affinity, Stoichiometry, Enthalpy	Affinity, Stoichiometry, Enthalpy
Tryptophan fluorescence required	No	No
Measurements in complex bioliquids (serum, cell lysate)	Yes	Yes
Volume per measurement	$< 4\mu\text{l}$	$< 4\mu\text{l}$
Molecular weight range (Da)	$10^1 - 10^7$	$10^1 - 10^7$
Time for experiment & analysis	Minutes	Minutes
Immobilization required	No	No
Temperature controlled	22 – 45°C	22 – 45°C
Maintenance required	No	No


Detectors and Spectra

The excitation/emission wavelengths for the corresponding Monolith NT.115 Series detectors are illustrated below.

The light and dark colored boxes correspond to the excitation and emission wavelength spectra, respectively, of the different fluorescence detectors.

Protein–Protein Interaction

Thermophoretic analysis of a protein–protein interaction


This interaction was investigated in a purified system using a covalently attached fluorophore (A) and using fluorescent fusion proteins in buffer (B) as well as in pure eukaryotic cell lysate (C). Comparable binding affinities were determined for this interaction in all three setups.

Material was kindly provided by Prof. Gideon Schreiber, Weizmann Institute, Israel

Jerabek–Willemsen, M., Andre, T., Wanner, R., Roth, H. M., Duhr, S., Baaske, P., and Breitsprecher, D. (2014)

MicroScale Thermophoresis: Interaction analysis and beyond. Journal of Molecular Structure

High-Affinity Interactions

Detection of picomolar K_d s of different experimental systems using the Monolith NT.115^{Pico}.

DNA–DNA interaction (left). The hybridization of a 16mer with a complementary Cy5–labeled DNA strand was monitored with MST. The change in MST signals was fitted (redline) to yield a K_d of $167 \pm 8 \text{ pM}$.

Protein–DNA interaction (center). Association of single strand binding protein from *E. coli* (EcoSSB) with a Cy5–labeled oligo–nucleotide (dT)70 was monitored by titrating EcoSSB against 20 pM Cy5–labeled DNA. Changes in thermophoresis were plotted, yielding a K_d of $1 \pm 0.1 \text{ pM}$.

Protein–protein interaction (right). Native Anakinra was titrated against NT647–labeled, soluble fragment of the Interleukin–1 receptor. Plotting of the change in thermophoresis and concomitant fitting of the data yielded a K_d of $620 \pm 61 \text{ pM}$.

Material was kindly provided by Dr. Ute Curth, Medical University Hanover, and Dr. Ahmed Besheer, Novartis, Basel.

Jerabek–Willemsen, M., Andre, T., Wanner, R., Roth, H. M., Duhr, S., Baaske, P., and Breitsprecher, D. (2014) MicroScale Thermophoresis: Interaction analysis and beyond. *Journal of Molecular Structure*

Monolith NT_LabelFree

Product Information

Monolith® Instruments
for MicroScale Thermophoresis

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
Monolith Instruments NT.Series	NT.115
Samples per run	16 samples
Fluorescence channels per instrument	1 (UV)
Fluorescence multiplexing	No
Affinity range	10 nM to mM
Labeling required	No
Concentration of fluorescent molecule	$10^{-8} - 10^{-3}$ M
Sample consumption	1.2 μ g ¹
Biophysical parameters	Affinity, Stoichiometry, Enthalpy
Tryptophan fluorescence required	Yes
Measurements in complex bioliquids (serum, cell lysate)	< 4 μ l
Volume per measurement	No
Molecular weight range (Da)	$10^1 - 10^7$ Da
Time for experiment & analysis	Minutes
Immobilization required	No
Temperature controlled	22 – 45°C
Maintenance required	No

¹calculated for a standard protein of 50 kDa, 12 data points per Kd and 100 nM tryptophan-containing protein.


Protein – Small Molecule Interaction

MST binding experiment of a GPCR Receptor binding to an inhibitor in the presence and absence of the natural ligand Neurotensin.

Prof. Anthony Watts, University of Oxford, Biochemistry, UK

Seidel et al, Angewandte Chemie, 2012


MST binding experiment of Synaptotagmin, a calcium sensor, to divalent cations.

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
² Max Planck Institute for Biophysical Chemistry, Göttingen, Germany

Protein – Sugar Interaction

MST experiment of maltose binding protein (MBP) binding to the sugars maltose and N, N', N'' – triacetylchitotriose.

MBP was kindly provided by Susanna v. Gronau and Dr. Sabine Suppmann from the biochemistry core facility of the Max Planck Institute of Biochemistry, Martinsried.

Interaction of Lytic Amidase LytA with two different mucopeptides

Mellroth, P., Sandalova, T., Kikhney, A., Vilaplana, F., Hesek, D., Lee, M., Mobashery, S., Normark, S., Svergun, D., Henriques-Normark, B., and Achour, A. (2014) Structural and Functional Insights into Peptidoglycan Access for the Lytic Amidase LytA of *Streptococcus pneumoniae*. *MBio* 5(1); doi:10.1128/mBio.01120-13.

amended with permission from the American Society for Microbiology

Prometheus[®] NT.48

Product Information

Prometheus Instruments[®]
for nanoDSF[®]

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Enjoy the benefits of nanoDSF :

- ▶ Measure thermal and chemical stability – even for membrane proteins
- ▶ See more transitions – ultra-high resolution for antibody engineering
- ▶ Detect aggregates – long-term stability and storage of biologics
- ▶ Exploit the concentration range – for formulation of biopharmaceuticals
- ▶ Integrate into robotic platforms – for fully automated operation

Prometheus NT.48


NanoDSF®

A technology by NanoTemper®

nanoDSF 는 advanced Differential Scanning Fluorimetry technology 입니다.

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
Thermal unfolding with nanoDSF

Two wavelengths, 330 and 350 nm are recorded. The ratio of the two wavelengths is plotted against the temperature. The 350/330 nm ratio typically yields well-defined transitions, even if the single wavelengths do not exhibit a clear unfolding transition.

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

Chemical Unfolding in Antibody Development

The Prometheus software package detects chemical unfolding and determines ΔG within seconds. Since concentration-dependence of ΔG is a measure for aggregation propensity, nanoDSF is able to predict long-term stability of biopharmaceuticals in different formulations.

The fully automated nanoDSF solution with the Prometheus NT,Plex allows for unattended measurement of hundreds of chemical unfolding reactions per day.


Determine ΔG of folding to predict aggregation propensities

Upper panel: Schematic representation of chemical unfolding experiments. Denaturant dilutions in formulation buffer are prepared in 384-well MTPs, and loaded into Prometheus capillary chips. Chemical denaturation curves are recorded and analyzed automatically in seconds.

Lower panel: Concentration-dependence of ΔG of a mAb in different formulations. A decrease in ΔG indicates high aggregation propensity of the unfolded state. Formulation 4 shows a constant ΔG , and also shows the highest monomer content after 12 months at 25°C by HPSEC.

Comparison: nanoDSF and μ DSC


A collaborative study with a leading CRO demonstrates that nanoDSF overcomes critical drawbacks associated with μ DSC since it provides perfect ease of use, is 100 times faster and requires 40 times less sample.

Both nanoDSF and μ DSC generate precise and highly comparable T_m values (Figure A and B) in a small formulation screen using a commercially available therapeutic mAb. A total of ten different formulations with varying buffers and pH-values were tested in addition to polysorbate 20 and 80, which are common surfactants in mAb formulations, but preclude the analysis by orthogonal fluorescence methods such as DSF assays.

The integrated nanoDSF protocol and innovative capillary sample format of the Prometheus NT.48 overcomes many key limitations of μ DSC (Figure C and D). In addition to its speed, precision and throughput, nanoDSF is a robust method that does not require cumbersome instrument maintenance and time-consuming sample preparation such as dialysis or filtration. Therefore, the Prometheus NT.48 is the ideal instrument for rapid and precise thermal stability screening in biopharmaceutical development.

Thermal and Colloidal Stability of Antibodies

Conformational stability and aggregation of a monoclonal antibody (mAb) under different buffer conditions.


(A) Thermal unfolding monitored by detection of shifts in the fluorescence ratio (F_{350}/F_{330}) in dependence of different buffer pH values and NaCl concentrations.

(B) Aggregation detected by changes in backreflection.

The Prometheus NT.48 delivers high quality thermal unfolding data for antibody buffer screening campaigns and detects unfolding of single antibody domains (A). The large dynamic range of the Prometheus NT.48 allows for analyzing thermal unfolding in solutions containing antibody concentrations between 250 mg/ml down to few μ g/ml. In addition, the colloidal stability (B) and aggregation onset temperatures (C) of antibodies can be assessed simultaneously employing the backreflection optics.

The results from our parallel investigation of thermal unfolding and aggregation using the Prometheus NT.48 suggest that a slight thermal instability might be acceptable and even favourable for long-term stability of the antibody, due to reduced aggregation under these conditions (D). Future screening approaches could therefore be designed to find excipients which thermally stabilize the mAb at lower pH values while maintaining low aggregation of the unfolded state.

Quality Control

Establishing a protein unfolding standard.

Unfolded IgG at different concentrations was mixed with folded IgG and subjected to thermal unfolding. The percentage of unfolded IgG in the solution was quantified based on the F350/F330 ratio measured at 25°C.

Forced-degradation stress-test on MEK1

MEK1 protein was subjected to the indicated stresses, and the fraction of unfolded protein was calculated based on the F350/F330 ratio at 25°C. Error bars are s.d. from three measurements.

nanoDSF can be employed to quickly detect and quantify unfolded proteins for quality control purposes with unmatched speed, at the same time offering unique ease of use.

The presented quality control experiments can be performed by filling capillaries directly from stock solutions without laborious sample preparation. F350/F330 values of up to 48 samples are then recorded in parallel using a one-button routine, providing stability data within seconds.

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